

## AED INSTRUCTION

### ZOLL AED PLUS

1. REMOVE THE GREEN COVER AND START THE DEFIBRILLATOR AT THE PURPLE ON/OFF BUTTON RECHTS BELOW.
2. PLACE THE ELECTRODES ACCORDING TO THE PICTURES ON THE NAKED AND DRY CHEST. IF MORE HELPERS ARE PRESENT, CONTINUE THE HEART PRESSURE MASSAGE WHILE PLACING THE ELECTRODES.
3. FOLLOW THE INSTRUCTIONS OF THE DEFIBRILLATOR. IF A SHOCK IS REQUIRED, PRESS THE BUTTON WITH THE LIGHTNING HEART SYMBOL IN THE MIDDLE OF THE DEFIBRILLATOR.

